



The medical lockable line valve comprises of a 2-piece full bore male threaded nickel plated brass ball valve c/w chrome plated brass ball, blow out proof stem, Viton® Stem O-ring, Teflon® ball seals, 2x Locknuts with visual indicator for blanking facility and flat face seal incorporating grade 316 Stainless steel stub pipe assemblies.

If In Doubt Ask							
Do Not Scale	Drawn by K.J.P.	Checked by	Approved by - date	Filename	Date 02/10/2014	Scale NTS	All Dims In mm Unless Otherwise Stated
	PRECISION UK Ltd			Lockable Line Valve Assembly ST/ST Stub Pipes			
	<small>CPX Technology Building, Pepper Road, Hazel Grove, Sockport SK7 5BW Tel: 00 44 (0) 161 487 2822 Fax: 00 44 (0) 161 487 2816 Email: info@precision-parts.co.uk www.precision-parts.co.uk</small>			LLV-ASS-3D-03-SS	Revision 0	Sheet 1/1	