

---

## MEDICAL PENDANTS - FIXED

---

### DESCRIPTION

Rigid and Retractable Pendants are manufactured under BS EN 13485 Medical Devices: Quality Management Systems. They are intended to distribute medical gases, AGSS and electrical services to the operating table/anaesthetic machine, whilst leaving the surrounding area free from trailing hoses and cables.

Rigid Pendants are ceiling mounted fixed columns complete with first fix plate and ceiling shroud, with an epoxy powder coated finish.

Retractable Pendants are ceiling mounted, height adjustable columns complete with first fix plate and ceiling shroud, with an epoxy powder coated finish. The height to the underside of the pendant can be adjusted throughout a 300mm pre-determined range.

A maximum of 9 medical gas and AGSS terminal units, complete with hose assemblies, can be accommodated within each pendant.

Each pendant can accommodate a maximum of 8 electrical sockets. Other specialist outlets such as data sockets, equipotential sockets, etc. can be incorporated into the pendant as required.



### DESIGNED TO MEET THE REQUIREMENTS OF

- HTM 02-01: Health Technical Memorandum Medical Gas Pipeline Systems
- BS EN 7396: Medical Gas Pipelines for Compressed Medical Gases and Vacuum
- BS EN ISO 11197: Medical Supply Units
- BS EN 15908: Non-interchangeable Screw Threaded (NIST) Low-pressure Connectors for Medical Gases
- BS EN 9170: Terminal Units for use with Compressed Medical Gases and Vacuum



# MEDICAL PENDANTS - FIXED

## ORDERING INFORMATION

Please state Rigid or Retractable Pendant and number of gases required. Use code RG for Rigid Pendant and RT for Retractable Pendant

**Example:**

PEND-RG-2 = Rigid Pendant – 2 Gas

PEND-RT-5 = Retractable Pendant – 5 Gas

