


RIGID PENDANT SPECIFICATION SHEET

Sales Contract No		Customer Order No.	
Project Ref.		Delivery Date Requested	
Qty		Batch No.	
Customer			

Gas Outlets	
Position	Gas Type
1	
2	
3	
4	
5	
6	
7	
8	
9	

Electrical Sockets	
Double Qty	
Switched or Unswitched	
Colour- White or Blue	
UPS circuit	
Ceiling Height	
Pendant Length	
Signed By	
Date	
Customer Signed	
Date	



It is stated that pendants should be clear 2000mm from floor to bottom of pendant e.g. 3mtr floor to ceiling height = 1mtr long pendant

Please return completed and signed off to enable production to begin, any delay may affect the agreed delivery date.

RIGID1/11