

Precision UK Ltd.
CPX Technology Building
Pepper Road
Hazel Grove
Cheshire
Sk7 5BW

Declaration of Conformity to EU RoHS

Products listed below that are manufactured by Precision UK Ltd are in compliance with Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment (also known as “RoHS Recast”). In addition, this declaration of conformity is issued under the sole responsibility of Precision UK Ltd. Specifically, products manufactured do not contain the substances listed in the table below in concentrations greater than the listed maximum value.

Substance	Maximum Limit (ppm)
Lead (Pb)	1000
Cadmium (Cd)	100
Mercury (Hg)	1000
Hexavalent Chromium (Cr6+)	1000
Poly Brominated Biphenyls (PBB)	1000
Poly Brominated Diphenyl ethers (PBDE)	1000

Product description:

The medical gas terminal unit comprises of a Second Fix assembly and Check Valve assembly which are fitted to a First Fix assembly. The first fix assembly includes a 12mm copper pipe for brazing to the fixed pipeline system or a non-interchangeable screw thread (NIST) connector for connection to a hose assembly.

www.precision-parts.co.uk



Company Registered Office: Crown House, Stockport, Cheshire, SK1 3RB
Registered in England and Wales No. 05058855
VAT No. 850 0700 67

Unit 8, Newby Rd Ind Est,
Stockport, Cheshire, UK, SK7 5DA
Tel: +44 (0) 161 487 2822
Fax: +44 (0) 161 487 2816
Email: info@precision-parts.co.uk

Product Identification:

Terminal Units - Flush Fitting	
Product Code	Gas Type
TUFO2	Oxygen
TUFN2O	Nitrous Oxide
TUFO2N2O	Entonox
TUFMA4	Medical Air 4
TUFMA7	Medical Air 7
TUFVAC	Vacuum
TUFAGSS	AGSS
TUFCO2	Carbon Dioxide
Terminal Units - Bedhead Fitting	
Product Code	Gas Type
TUBHO2	Oxygen
TUBHN2O	Nitrous Oxide
TUBHO2N2O	Entonox
TUBHMA4	Medical Air 4
TUBHMA7	Medical Air 7
TUBHVAC	Vacuum
TUBHAGSS	AGSS
Terminal Units - Back Entry NIST	
Product Code	Gas Type
TUBEO2	Oxygen
TUBEN2O	Nitrous Oxide
TUBEO2N2O	Entonox
TUBEMA4	Medical Air 4
TUBEMA7	Medical Air 7
TUBEVAC	Vacuum
TUBEAGSS	AGSS

Date of Issue: 20/10/2015

Signature:

Name (printed): M Dean

Title: Engineer

www.precision-parts.co.uk



Company Registered Office: Crown House, Stockport, Cheshire, SK1 3RB
Registered in England and Wales No. 05058855
VAT No. 850 0700 67

Unit 8, Newby Rd Ind Est,
Stockport, Cheshire, UK, SK7 5DA
Tel: +44 (0) 161 487 2822
Fax: +44 (0) 161 487 2816
Email: info@precision-parts.co.uk